

SPH Series of Hydrogen Generators



Model	SPH-200	SPH-300A	SPH-300	SPH-500A	SPH-500
Specifications					
Hydrogen Purity	99.999%	99.999%	99.999%	99.999%	99.999%
Hydrogen Flow	0-200ml/min	0-300ml/min	0-300ml/min	0-500ml/min	0-500ml/min
Working Pressure	0.3MPa	0.4MPa	0.4MPa/0.6MPa (2-step)	0.4MPa	0.4MPa/0.6MPa (2-step)
Power Consumption	100W	150W	150W	250W	250W
Dimension (mm)	200x140x290	360x200x260	370x180x330	360x200x260	370x180x330
Net Weight	6kg	10kg	10kg	10kg	10kg

Features:

1. Fully automatic
2. No backward alkali solution